



Max Force partial foot prosthesis User instructions:

The Max Force system has allowed your clinical prosthetist to create today's most advanced partial foot prosthesis. It is custom made to the mold of your lower leg and built to the exact specifications of your clinician. Your clinician is critical to achieving your maximum performance with the Max Force system. You should always contact your prosthetist first if you have any questions or concerns about the fit or function of your Max Force Prosthesis.

Always !

Wear a sock on your lower leg. For the most comfort, the sock should extend above the prosthetic socket and should not have large wrinkles.

Always !

Wear a shoe when using the Max Force Prosthesis. You will likely slip and fall if you walk without a proper shoe.

Always !

Inspect your leg if you feel discomfort. If this continues or if you develop any blisters or abrasions of the skin, contact your prosthetic clinician.

Always !

Inspect your leg regularly to ensure that you have not developed a blister or lesion that you may not be aware of. This is especially important if you do not have normal sensation in your leg

Cleaning and Maintenance

**Hand wash with mild detergent and water throughout the device as needed.
air dry only**

No specific maintenance is required. however, it is recommended that pads and straps be replaced on a regular basis, as needed.