



# Instructions for T-Strap Use



- 
- 1) You will receive your T-Strap with a small strip of sticky back hook Velcro.
- 

- 
- 2) Place the strip of Velcro on the distal tab of the T-Strap then remove the plastic backing.
- 



- 
- 3) Align the T-Strap at the appropriate placement for your patient and attach the Velcro to the bottom of the footplate.
- 



- 
- 4) Place the Noodle and T-Strap into the patient's shoe.
- 

- 5) Have the patient don the AFO and fasten the proximal strap.
- 

- 6) Loop the ankle strap around the strut and through the D-Ring. Tighten to get varus/valgus correction.

