



Instructions for T-Strap Use



-
- 1) You will receive your T-Strap with a small strip of sticky back hook Velcro.
-

-
- 2) Place the strip of Velcro on the distal tab of the T-Strap then remove the plastic backing.
-



-
- 3) Align the T-Strap at the appropriate placement for your patient and attach the Velcro to the bottom of the footplate.
-



-
- 4) Place the Noodle and T-Strap into the patient's shoe.
 - 5) Have the patient don the AFO and fasten the proximal strap.
 - 6) Loop the ankle strap around the strut and through the D-Ring. Tighten to get varus/valgus correction.
-

