



# Custom From Cast Noodle AFO Order Form

5513 WEST SLIGH AVE. TAMPA, FL 33634

PHONE: 800-919-3668 FAX: 888-484-5927 WWW.KINETICRESEARCH.COM EMAIL: EMAIL@KINETICR.COM

Company: \_\_\_\_\_

Location: \_\_\_\_\_

PO#: \_\_\_\_\_

Date: \_\_\_\_\_ Date due: \_\_\_\_\_

Practitioner: \_\_\_\_\_

Patient name: \_\_\_\_\_

Age: \_\_\_\_\_ Height: \_\_\_\_\_ Weight: \_\_\_\_\_

Male  Female

Activity level: \_\_\_\_\_

Pathology: \_\_\_\_\_

Cast position \_\_\_\_\_

Correct to \_\_\_\_\_

All casts will be corrected to 90 degrees on a 3/8" heel height and subtalar neutral unless otherwise specified.

**Left**  **Right**

### Noodle AFO Style:

Cuff: Posterior Classic  
Posterior Hemi-spiral  
Anterior TA

### Footplate:

Flexibility:

Flexible

**Normal**

Stiff

Extra stiff

Style:

**Flat:**



**Contoured\*:**  
(must order with extra stiff strut)



**Strut:** Medial

**Lateral**

Flexible

**Normal**

Stiff

Extra Stiff

Heavy Duty Layup\*

### Practitioner Notes and Drawspace

AFO Height: \_\_\_\_\_

Footplate Length: \_\_\_\_\_

\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_

### Shipping:

- Ground
- 3-day
- 2 day
- Next day PM
- Next Day AM

### Pre-Tibial Shell\*:

(for use with posterior cuff)

### Long Anterior Cuff\*:

#### Molded Foot Orthotic\*:

Diabetic Insert  
Cork/EVA Rough  
Cork

#### Molded Inner Boot\*:

UCBL: 1/8" Polypropylene  
SMO: 1/8" Polyethylene

### T-Strap\*:

### Color:

Black

Caucasian\*

Fabric\*: \_\_\_\_\_